## DEUTSCH

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## ENGLISH

### Important safety information – retain for later use.

In order to make optimum use of all the merits of the clinical thermometer, you should carefully read the operating instructions before using it, retain the instructions for later use and also keep them accessible for other users.

- The thermometer is intended to measure the temperature of the human body.
- Children should not use the thermometer without supervision.
- Before each use, check the thermometer for signs of damage or wear. Damaged or worn thermometers should not be used.
- Always wait for the acoustic signal which indicates the minimum time required for measurement. Discuss the temperature with your family physician.

### Speicher

The last measured value will be automatically stored. This value will be displayed when the EIN/AUS button is pressed. In addition to the measured value, the storage symbol, e.g. 37,5°C (99,5°F), will appear.

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You can change the temperature unit to Grad Celsius (°C) and Grad Fahrenheit (°F) by displaying the value. Hold the EIN/AUS button until the temperature value is displayed for longer than 2 seconds. The measured value will then be displayed in the new unit.

### Batteriewechsel

The clinical thermometer contains a long-life alkali-manganese battery (LR 41, 1,55 V), which will last for approx. 2 years. The battery has to be replaced when 'T' appears in the bottom right corner of the display. Press down the battery cover at the end of the device with your fingers and then pull it firmly back. Remove the old battery and replace it with a battery of the same type. The + sign on the battery must point upwards. Then put the battery cover back onto the case. Used batteries do not belong in the regular domestic refuse. Dispose of them through your local dealer for electric goods or your local depot for special waste. You are legally obliged to dispose of the batteries. Note: Batteries containing pollutant substances are marked as follows:

Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

### Cleaning/Disinfection

The thermometer is 100% water-proof. For the purpose of cleaning or disinfection, it can be dipped into water or disinfectant without any difficulties. Under no circumstances should the thermometer be boiled! Use one of the following disinfectants in accordance with the manufacturer's instructions. Other cleaning agents or methods could impair the function of the thermometer or damage it.

**Name of the agent:** Isopropanol 70% Pharmacy; Gigasept FF; Lysoformin; Aseptisol.

### Storage/Disposal

When you are not using the thermometer, please retain it in the original protective cover. The device must be tested for accuracy (calibrated) after 2 years. Please send the complete device to the dealer's or service address. Please dispose of the device in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

### Technical data

Type: Maximum thermometer type SFT 11/1

Range of measurement: 32 °C (89,6°F) to 42,9 °C (109,2°F)  
Accuracy of measurement: ±0,1 °C in a water bath with a temperature of 35,5 °C (95,9°F) to 42,0 °C (107,8°F)

±0,2 °C in a water bath with a temperature of 32,0 °C (89,6°F) to 35,4 °C (95,7°F)  
±0,2 °C in a water bath with a temperature of 42,0 °C (107,8°F) to 42,9 °C (109,2°F)

Ambient temperature: 10 °C (50°F) to 40 °C (122°F), with 30% to 85% relative humidity

Storage temperature: -10 °C (14°F) to +60 °C (140°F), with 25% to 90% relative humidity

Reference of standardisation: EN 12470-3: 2000 Clinical Thermometers

Explanation of symbols: Attention – read the operating instructions. Application part of Type BF

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